

## FIRST DRAW & FLUSH LEAD TEST SAMPLING INSTRUCTIONS

**PLEASE READ ALL INSTRUCTIONS COMPLETELY BEFORE COLLECTING YOUR WATER SAMPLE**

**One test kit represents one sampling location. If you require testing on multiple locations you will need to purchase multiple test kits.**

1. First Draw sample: Choose a collection time when the water has not been used for at least six (6) hours. Before you use any water from the tap, fill the sample vial labeled "First Draw" as soon as you turn on the tap.
2. Flush sample: Run the water from the tap for 3-5 minutes and then fill the sample vial labeled "Flush".
3. Complete the lab form, ***making sure to indicate the date and time you sampled***, type of water sampled (treated/untreated) and sample identification which identifies where you collected your sample (i.e., kitchen sink, well head, etc.). Please be sure to retain a copy of the completed lab form prior to shipment. Return the completed paperwork and water samples to the shipping container.
4. Mail your water sample to the testing location below. The sample does not require overnight shipping.

Testing Location: **National Testing Laboratories, Ltd.  
556 South Mansfield Street  
Ypsilanti, MI 48197**

5. PLEASE NOTE: Testing takes approximately 5-7 business days from receipt of the samples. Results can be emailed or faxed (if provided an email address or fax number).

If you have any questions regarding this First Draw & Flush Lead Test or any other services we provide, please contact a Technical Service Representative at National Testing Laboratories at 1-800-458-3330 for assistance.

Thank you for choosing National Testing Laboratories, Ltd., for your water testing services!