

QUICK SCREEN TEST SAMPLING INSTRUCTIONS

***Quick Screen is a Quick Series testing package designed to provide rapid test results.
This testing package includes Alkalinity, Iron, Hardness, Total Dissolved Solids and pH.
Results in 2-3 business days***

PLEASE READ ALL INSTRUCTIONS COMPLETELY BEFORE COLLECTING YOUR WATER SAMPLE

One test kit represents one sampling location. If you require testing on multiple locations you will need to purchase multiple test kits.

1. Freeze the blue gel ice pack for at least 8 hours prior to sample collection. Return to the cooler just prior to sampling and/or return shipping.
2. Remove aerators, strainers, mixers or purifiers from your tap.
3. Run the cold water for approximately 3-5 minutes, then adjust the flow to a thin stream to avoid splashing or overfilling.
4. Fill the large round bottle to the neck, cap securely and return it to the cooler.
5. Complete the enclosed lab form, **making sure to indicate the date and time you sampled**, type of water (treated/untreated) and sample identification which identifies where you collected your sample (i.e., kitchen sink, well head, etc.). Retain a copy of the lab form for your records and return the completed form to the shipping container.
6. Mail your sample to the testing location below. The sample should be sent via 2 - 3 day delivery.

Testing Location: **National Testing Laboratories, Ltd.
556 South Mansfield Street
Ypsilanti, MI 48197**

7. PLEASE NOTE: Testing takes approximately 2-3 business days from receipt of the samples. Results can be emailed or faxed (if provided an email address or fax number).

If you have any questions regarding this Quick Screen Test or any other services we provide, please contact a Technical Service Representative at National Testing Laboratories at 1-800-458-3330 for assistance.

Thank you for choosing National Testing Laboratories, Ltd., for your water testing services!